

Le changement passe par vous
Be part of the solution

FILLWARE PRINTING INSTRUCTIONS

- 1- To print an Ecopill label (compliance package), you go into a patient's file and click on **REPORTS**.
- 2- Scroll down and choose **ECOLOPHARM ECOPILL REPORT**.
- 3- You can then customize the look of the label to your liking. For example, if you want meds to be organized in a grid on the header, simply click on the grid option etc.

We suggest you play around with your options on plain paper to get an idea of the template you wish to work with. You have 6-7 options with a 14 med maximum.

How to place labels in printer tray:

Always place paper upside down (liner facing up) with the header towards you in the printer.

For double sided printing:

Print label as described above then replace same printed label with the header towards you but this time the coloured and printed part must be facing upwards.